

Antibodies Product Datasheet

Product Name:	AMACR Monoclonal Antibody(4A12)
Catalog Number:	N1311
Reactivity:	Human,Mouse,Rat
Applications:	WB,IF/ICC,IHC-p
Dilutions:	WB: 1:1000 IHC: 1:200 IF 1:200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of AMACR
Concentration:	
Gene Name:	AMACR
Protein Name:	Alpha-methylacyl-CoA racemase
Other Name:	AMACR; Alpha-methylacyl-CoA racemase; 2-methylacyl-CoA racemase
Human Swiss Prot No.:	Q9UHK6
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous AMCAR/P504S proteins.
Storage and Stability:	-20°C/1 year