

Antibodies Product Datasheet

Product Name:	CD68 Monoclonal Antibody(6F3)
Catalog Number:	N1314
Reactivity:	Human
Applications:	IHC-P,IF/ICC
Dilutions:	WB 500-2000 1:200 IF 1:50-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of CD68
Concentration:	
Gene Name:	CD68
Protein Name:	Macrosialin
Other Name:	CD68; Macrosialin; Gp110; CD68
Human Swiss Prot No.:	P34810
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous CD68 proteins.
Storage and Stability:	-20°C/1 year