

Antibodies Product Datasheet

Product Name:	CK7 Monoclonal Antibody(12D7)
Catalog Number:	N1319
Reactivity:	Human,Mouse,Rat
Applications:	IF/ICC,WB,IHC-pIP
Dilutions:	IF: 1:50-200 WB: 1:1000-2000 IHC: 1:200 IP:1:200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of CK7
Concentration:	
Gene Name:	KRT7
Protein Name:	Keratin type II cytoskeletal 7
Other Name:	KRT7; SCL; Keratin, type II cytoskeletal 7; Cytokeratin-7; CK-7; Keratin-7; K7; Sarcolectin; Type-II keratin Kb7
Human Swiss Prot No.:	P08729
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous CK7 proteins.
Storage and Stability:	-20°C/1 year