

Antibodies Product Datasheet

Product Name:	GFAP Monoclonal Antibody(5C8)
Catalog Number:	N1324
Reactivity:	Human,Rat,Mouse
Applications:	WB,IHC-p,IF/ICC
Dilutions:	WB: 1:2000-5000 IF 1:200 IHC 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of GFAP
Concentration:	
Gene Name:	GFAP
Protein Name:	Glial fibrillary acidic protein
Other Name:	GFAP; Glial fibrillary acidic protein; GFAP
Human Swiss Prot No.:	P14136
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous GFAP proteins.
Storage and Stability:	-20°C/1 year