

Antibodies Product Datasheet

Product Name:	CD21 Monoclonal Antibody(2C5)
Catalog Number:	N1325
Reactivity:	Human, Mouse, Rat
Applications:	IHC-P, IF/ICC
Dilutions:	IHC 1:200 IF 1:50-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of CD21
Concentration:	
Gene Name:	CR2
Protein Name:	Complement receptor type 2
Other Name:	CR2; C3DR; Complement receptor type 2; Cr2; Complement C3d receptor; Epstein-Barr virus receptor; EBV receptor; CD21
Human Swiss Prot No.:	P20023
Formulation:	PBS, pH 7.4, containing 0.5% BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous CD21 proteins.
Storage and Stability:	-20°C/1 year