

Antibodies Product Datasheet

Product Name:	Ki 67 Monoclonal Antibody(4A8)
Catalog Number:	N1329
Reactivity:	Human,Mouse,Rat
Applications:	IHC-P,IF/ICC
Dilutions:	IHC 1:200 IF 1:50-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of Ki 67
Concentration:	
Gene Name:	MKI67
Protein Name:	Antigen KI-67
Other Name:	MKI67; Antigen KI-67
Human Swiss Prot No.:	P46013
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous Ki 67 proteins.
Storage and Stability:	-20°C/1 year