

Antibodies Product Datasheet

Product Name:	MAP2 Monoclonal Antibody(7D4)
Catalog Number:	N1332
Reactivity:	Human,Mouse,Rat
Applications:	IHC-P,IF/ICC
Dilutions:	WB 500-2000 1:200 IF 1:50-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of MAP2
Concentration:	
Gene Name:	MAP2
Protein Name:	Microtubule-associated protein 2
Other Name:	MAP2; Microtubule-associated protein 2; MAP-2
Human Swiss Prot No.:	P11137
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous MAP2 proteins.
Storage and Stability:	-20°C/1 year