

Antibodies Product Datasheet

Product Name:	CD1 Monoclonal Antibody(9H6)
Catalog Number:	N1333
Reactivity:	Human,Mouse,Rat
Applications:	IHC-P,IF/ICC
Dilutions:	WB 500-2000 1:200 IF 1:50-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of CD1
Concentration:	
Gene Name:	CD1A
Protein Name:	T-cell surface glycoprotein CD1a
Other Name:	T-cell surface glycoprotein CD1a (T-cell surface antigen T6/Leu-6) (hTa1 thymocyte antigen) (CD antigen CD1a)
Human Swiss Prot No.:	P06126
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous CD1 proteins.
Storage and Stability:	-20°C/1 year