

Antibodies Product Datasheet

Product Name:	CD10 Monoclonal Antibody(5B8)
Catalog Number:	N1337
Reactivity:	Human,Mouse,Rat
Applications:	IF/ICC,IHC-p
Dilutions:	IF: 1:50-200 WB 500-2000 : 1:200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of CD10
Concentration:	
Gene Name:	MME
Protein Name:	Neprilysin
Other Name:	MME; EPN; Neprilysin; Atriopeptidase; Common acute lymphocytic leukemia antigen; CALLA; Enkephalinase; Neutral endopeptidase 24.11; NEP; Neutral endopeptidase; Skin fibroblast elastase; SFE; CD10
Human Swiss Prot No.:	P08473
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous CD10 proteins.
Storage and Stability:	-20°C/1 year