

Antibodies Product Datasheet

Product Name:	FH Monoclonal Antibody(7F1)
Catalog Number:	N1338
Reactivity:	Human,Mouse,Rat
Applications:	WB,IHC-p,IF/ICC
Dilutions:	WB: 1:3000 IF 1:200 IHC 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of FH
Concentration:	
Gene Name:	FH
Protein Name:	Fumarate hydratase, mitochondrial
Other Name:	Fumarate hydratase, mitochondrial (Fumarase) (EC 4.2.1.2)
Human Swiss Prot No.:	P07954
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous FH proteins.
Storage and Stability:	-20°C/1 year