

Antibodies Product Datasheet

Product Name:	AFP alpha 1 Fetoprotein Monoclonal Antibody(17C5)
Catalog Number:	N1340
Reactivity:	Human
Applications:	WB,IHC-p,IF/ICC
Dilutions:	WB: 1:2000 IF 1:200 IHC 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of AFP alpha 1 Fetoprotein
Concentration:	
Gene Name:	AFP
Protein Name:	Alpha-fetoprotein
Other Name:	AFP; HPAFP; Alpha-fetoprotein; Alpha-1-fetoprotein; Alpha-fetoglobulin
Human Swiss Prot No.:	P02771
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous AFP proteins.
Storage and Stability:	-20°C/1 year