

Antibodies Product Datasheet

Product Name:	Caspase 9 Monoclonal Antibody(3-20)
Catalog Number:	N1342
Reactivity:	Human,Mouse,Rat,chicken
Applications:	WB,IHC-p,IF/ICCIP
Dilutions:	WB: 1:1000-5000 IP:1:200 IF 1:200 IHC 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of Caspase 9
Concentration:	
Gene Name:	CASP9
Protein Name:	Caspase9
Other Name:	CASP9; MCH6; Caspase-9; CASP-9; Apoptotic protease Mch-6; Apoptotic protease-activating factor 3; APAF-3; ICE-like apoptotic protease 6; ICE-LAP6
Human Swiss Prot No.:	P55211
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous Caspase 9 protein.
Storage and Stability:	-20°C/1 year