

Antibodies Product Datasheet

Product Name:	ERCC1 Monoclonal Antibody(1B10)
Catalog Number:	N1343
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),
Dilutions:	IHC-p: 100-300.WB: 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of ERCC1
Concentration:	
Gene Name:	ERCC1
Protein Name:	DNA excision repair protein ERCC-1
Other Name:	ERCC1; DNA excision repair protein ERCC-1
Human Swiss Prot No.:	P07992
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous ERCC1 proteins.
Storage and Stability:	-20°C/1 year