

Antibodies Product Datasheet

Product Name:	IDE Monoclonal Antibody(3H4)
Catalog Number:	N1348
Reactivity:	Human
Applications:	WB,IHC-p,IF/ICC
Dilutions:	WB: 1:1000 IF 1:200 IHC 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of IDE
Concentration:	
Gene Name:	IDE
Protein Name:	Insulin-degrading enzyme
Other Name:	IDE; Insulin-degrading enzyme; Abeta-degrading protease; Insulin protease; Insulinase; Insulysin
Human Swiss Prot No.:	P14735
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous IDE proteins.
Storage and Stability:	-20°C/1 year