

Antibodies Product Datasheet

Product Name:	EFHD1 Monoclonal Antibody(3G2)
Catalog Number:	N1350
Reactivity:	Human,Mouse,Rat
Applications:	WB,IHC-p,IF/ICC
Dilutions:	WB: 1:2000 IF: 1:100-200 IHC 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of EFHD1
Concentration:	
Gene Name:	EFHD1
Protein Name:	EF-hand domain-containing protein D1
Other Name:	EF-hand domain-containing protein D1 (EF-hand domain-containing protein 1) (Swiprosin-2)
Human Swiss Prot No.:	Q9BUP0
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous EFHD1 proteins.
Storage and Stability:	-20°C/1 year