

Antibodies Product Datasheet

Product Name:	CD16 Monoclonal Antibody(Q32)
Catalog Number:	N1355
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),
Dilutions:	WB: 1:1000 IHC 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of CD16
Concentration:	
Gene Name:	FCGR3A/FCGR3B
Protein Name:	Low affinity immunoglobulin gamma Fc region receptor III-A/Low affinity immunoglobulin gamma Fc region receptor III-B
Other Name:	
Human Swiss Prot No.:	P08637/O75015
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous CD16 protein.
Storage and Stability:	-20°C/1 year