

Antibodies Product Datasheet

Product Name:	CD15 Monoclonal Antibody(Q89)
Catalog Number:	N1370
Reactivity:	Human
Applications:	IHC-P,IF/ICC
Dilutions:	IHC 1:200 IF 1:50-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of CD15
Concentration:	
Gene Name:	FUT4
Protein Name:	Alpha-(1,3)-fucosyltransferase
Other Name:	FUT4; ELFT; FCT3A; Alpha-(1,3)-fucosyltransferase; ELAM-1 ligand fucosyltransferase; Fucosyltransferase 4; Fucosyltransferase IV; Fuc-TIV; FucT-IV; Galactoside 3-L-fucosyltransferase
Human Swiss Prot No.:	P22083
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous CD15 protein.
Storage and Stability:	-20°C/1 year