

Antibodies Product Datasheet

Product Name:	CD41 Monoclonal Antibody(Q90)
Catalog Number:	N1371
Reactivity:	Human
Applications:	IHC-p,IF(paraffin section)
Dilutions:	IHC: 1:100
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of CD41
Concentration:	
Gene Name:	ITGA2B
Protein Name:	Integrin alpha-IIb
Other Name:	ITGA2B; GP2B; ITGAB; Integrin alpha-IIb; GPalpha IIb; GPIIb; Platelet membrane glycoprotein IIb; CD41
Human Swiss Prot No.:	P08514
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous human CD41/ Intergrin α 2b protein.
Storage and Stability:	-20°C/1 year