

Antibodies Product Datasheet

Product Name:	Peroxiredoxin 1 Monoclonal Antibody(8E7)
Catalog Number:	N1377
Reactivity:	Human,Mouse,Rat
Applications:	WB,IHC-p,IF/ICC
Dilutions:	WB: 1:1000-3000 IF: 1:100-200 IHC 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein of Peroxiredoxin-1
Concentration:	
Gene Name:	PRDX1
Protein Name:	Peroxiredoxin-1
Other Name:	PRDX1; PAGA; PAGB; TDPX2; Peroxiredoxin-1; Natural killer cell-enhancing factor A; NKEF-A; Proliferation-associated gene protein; PAG; Thioredoxin peroxidase 2; Thioredoxin-dependent peroxide reductase 2
Human Swiss Prot No.:	Q06830
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous Peroxiredoxin 1 protein.
Storage and Stability:	-20°C/1 year

N1377 | Peroxiredoxin 1 Monoclonal Antibody(8E7)

Tech@assaybiotech.com | Order@assaybiotech.com

Research use only not intended for therapeutic/diagnostic procedures

Last updated December 2022