

Antibodies Product Datasheet

Product Name:	Insulin Glargine Monoclonal Antibody(3F12)
Catalog Number:	N1381
Reactivity:	Species independent
Applications:	ELISA
Dilutions:	ELISA: 1:2000-10000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Insulin Glargine
Concentration:	
Gene Name:	
Protein Name:	
Other Name:	
Human Swiss Prot No.:	
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects Insulin Glargine protein.
Storage and Stability:	-20°C/1 year