

Antibodies Product Datasheet

Product Name:	Collagen III Monoclonal Antibody(Q76)
Catalog Number:	N1388
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC
Dilutions:	WB: 1:1000 IF 1:200 IHC 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of Collagen III
Concentration:	
Gene Name:	COL3A1
Protein Name:	Collagen alpha-1(III) chain
Other Name:	COL3A1; Collagen alpha-1(III) chain
Human Swiss Prot No.:	P02461
Formulation:	PBS, pH 7.4, containing 0.5% BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous Collagen III protein.
Storage and Stability:	-20°C/1 year