

Antibodies Product Datasheet

Product Name:	eIF4A1 Monoclonal Antibody(M8)
Catalog Number:	N1401
Reactivity:	Human,Mouse,Rat
Applications:	WB,IHC-p,IF/ICC
Dilutions:	WB: 1:1000-3000 IF: 1:100-200 IHC 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of eIF4A1
Concentration:	
Gene Name:	EIF4A1
Protein Name:	Eukaryotic initiation factor 4A-I
Other Name:	Eukaryotic initiation factor 4A-I (eIF-4A-I) (eIF4A-I) (EC 3.6.4.13) (ATP-dependent RNA helicase eIF4A-1)
Human Swiss Prot No.:	P60842
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous eIF4A1 protein.
Storage and Stability:	-20°C/1 year