

Antibodies Product Datasheet

Product Name:	Fibronectin Monoclonal Antibody(M9)
Catalog Number:	N1402
Reactivity:	Human,Mouse,Rat
Applications:	WB,IHC-p,IF/ICC
Dilutions:	WB: 1:1000-2000 IF 1:200 IHC 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of Fibronectin
Concentration:	
Gene Name:	FN1
Protein Name:	Fibronectin
Other Name:	FN1; FN; Fibronectin; FN; Cold-insoluble globulin; CIG
Human Swiss Prot No.:	P02751
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous Fibronectin protein.
Storage and Stability:	-20°C/1 year