

Antibodies Product Datasheet

Product Name:	Cleaved PARP Monoclonal Antibody(Mix)
Catalog Number:	N1410
Reactivity:	Human
Applications:	IF/ICC, WB, IHC-p
Dilutions:	IF: 1:50-200 WB: 1:2000-5000 IHC 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of Cleaved PARP
Concentration:	
Gene Name:	PARP1
Protein Name:	Poly [ADP-ribose] polymerase 1
Other Name:	PARP1; ADPRT; PPOL; Poly [ADP-ribose] polymerase 1; PARP-1; ADP-ribosyltransferase diphtheria toxin-like 1; ARTD1; NAD(+) ADP-ribosyltransferase 1; ADPRT 1; Poly[ADP-ribose] synthase 1
Human Swiss Prot No.:	P09874
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous pro and active PARP protein.
Storage and Stability:	-20°C/1 year