

Antibodies Product Datasheet

Product Name:	Transferrin Monoclonal Antibody(7F4)
Catalog Number:	N1417
Reactivity:	Human
Applications:	WB,IHC-p,IF/ICC
Dilutions:	WB: 1:1000-2000 IF 1:200 IHC 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of Transferrin
Concentration:	
Gene Name:	TF
Protein Name:	Serotransferrin
Other Name:	TF; Serotransferrin; Transferrin; Beta-1 metal-binding globulin; Siderophilin
Human Swiss Prot No.:	P02787
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous Human Transferrin protein.
Storage and Stability:	-20°C/1 year