

Antibodies Product Datasheet

Product Name:	Lactoferrin Monoclonal Antibody(Q100)
Catalog Number:	N1425
Reactivity:	Human
Applications:	IF/ICC,IHC-p,ELISA
Dilutions:	IF: 1:50-200 ELISA: 1:5000-10000 IHC 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of Lactoferrin
Concentration:	
Gene Name:	LTF
Protein Name:	Lactotransferrin
Other Name:	LTF; LF; Lactotransferrin; Lactoferrin; Talalactoferrin
Human Swiss Prot No.:	P02788
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous Lactoferrin protein.
Storage and Stability:	-20°C/1 year