

Antibodies Product Datasheet

Product Name:	mStrawberry Monoclonal Antibody(Mix)
Catalog Number:	N1428
Reactivity:	Species independent
Applications:	WB
Dilutions:	WB: 1:2000-5000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of mStrawberry
Concentration:	
Gene Name:	
Protein Name:	
Other Name:	
Human Swiss Prot No.:	
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects mStrawberry & mStrawberry Tag protein. It May had cross reactivity with mRFP protein.
Storage and Stability:	-20°C/1 year