

Antibodies Product Datasheet

Product Name:	eNOS Monoclonal Antibody(Mix)
Catalog Number:	N1429
Reactivity:	Human,Mouse,Rat
Applications:	WB
Dilutions:	WB: 1:500-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein of eNOS
Concentration:	
Gene Name:	NOS3
Protein Name:	Nitric oxide synthase, endothelial
Other Name:	NOS3; Nitric oxide synthase, endothelial; Constitutive NOS; cNOS; EC-NOS; Endothelial NOS; eNOS; NOS type III; NOSIII
Human Swiss Prot No.:	P29474
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous eNOS protein.
Storage and Stability:	-20°C/1 year