

Antibodies Product Datasheet

Product Name:	Rabbit IgG (Fc specific, Heavy Chain) Monoclonal Antibody(1E1)
Catalog Number:	N1432
Reactivity:	Rabbit
Applications:	WB
Dilutions:	WB: 1:2000-5000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein of Rabbit IgG (Fc specific, Heavy Chain)
Concentration:	
Gene Name:	
Protein Name:	
Other Name:	
Human Swiss Prot No.:	
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody reacts with the Fc portion of rabbit IgG heavy chain but not with the Fab portion of rabbit immunoglobulins. The antibody may cross-react with immunoglobulins from other species.
Storage and Stability:	-20°C/1 year