

Antibodies Product Datasheet

Product Name:	Rabbit IgG (Light Chain Specific) Monoclonal Antibody(8D6)
Catalog Number:	N1433
Reactivity:	Rabbit
Applications:	WB
Dilutions:	WB: 1:2000-5000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein of Rabbit IgG (Light Chain Specific)
Concentration:	
Gene Name:	
Protein Name:	
Other Name:	
Human Swiss Prot No.:	
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody reacts with the light chains on rabbit IgG and with those common to other rabbit immunoglobulins. The antibody does not react with the heavy chain of rabbit IgG. The antibody minimal cross-reaction with bovine, goat, Armenian hamster, horse,
Storage and Stability:	-20°C/1 year