

## Antibodies Product Datasheet

|                               |                                                                                                            |
|-------------------------------|------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | EYFP Monoclonal Antibody(Mix)                                                                              |
| <b>Catalog Number:</b>        | N1435                                                                                                      |
| <b>Reactivity:</b>            | Species independent                                                                                        |
| <b>Applications:</b>          | WB,IF                                                                                                      |
| <b>Dilutions:</b>             | WB: 1:5000,IF:1:300-500                                                                                    |
| <b>Modification:</b>          | Unmodified/ Total                                                                                          |
| <b>Source:</b>                | Monoclonal Mouse                                                                                           |
| <b>Purification:</b>          | The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen. |
| <b>Immunogen:</b>             | Recombinant Protein of EYFP                                                                                |
| <b>Concentration:</b>         |                                                                                                            |
| <b>Gene Name:</b>             |                                                                                                            |
| <b>Protein Name:</b>          |                                                                                                            |
| <b>Other Name:</b>            |                                                                                                            |
| <b>Human Swiss Prot No.:</b>  |                                                                                                            |
| <b>Formulation:</b>           | PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.                      |
| <b>Specificity:</b>           | The antibody detects EYFP and YFP tag fusion proteins.                                                     |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                               |