

Antibodies Product Datasheet

Product Name:	C. Elegans α -tubulin Mouse mAb(Mix-mA)
Catalog Number:	N1443
Reactivity:	C. Elegans (H, M, R)
Applications:	IHC-p,WB
Dilutions:	IHC-p1:100-500 ,WB 1:5000-20000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from C. Elegans α -tubulin
Concentration:	1 mg/ml
Gene Name:	C. Elegans α -tubulin
Protein Name:	C. Elegans α -tubulin
Other Name:	
Human Swiss Prot No.:	H2L2E8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.09% sodium azide.
Specificity:	This antibody detects endogenous levels of C. Elegans α -tubulin
Storage and Stability:	-20°C/1 year