

Antibodies Product Datasheet

Product Name:	TRP2 DCT Mouse mAb(2A10)
Catalog Number:	N1455
Reactivity:	Human, Mouse,Rat
Applications:	IHC-p,WB
Dilutions:	IHC-p1:50-200 ,WB 1:1000-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human TRP2 DCT
Concentration:	1 mg/ml
Gene Name:	DCT TYRP2
Protein Name:	TRP2 DCT
Other Name:	L-dopachrome tautomerase (DCT) (DT) (EC 5.3.3.12) (L-dopachrome Delta-isomerase) (Tyrosinase-related protein 2) (TRP-2) (TRP2)
Human Swiss Prot No.:	P40126
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.28% sodium azide.
Specificity:	This antibody detects endogenous levels of TRP2 DCT at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year