

Antibodies Product Datasheet

Product Name:	Calretinin Mouse mAb(Mix-mA)
Catalog Number:	N1456
Reactivity:	Human, Mouse,Rat
Applications:	IHC
Dilutions:	IHC 1:50-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Calretinin
Concentration:	1 mg/ml
Gene Name:	CALB2 CAB29
Protein Name:	Calretinin
Other Name:	Calretinin (CR) (29 kDa calbindin)
Human Swiss Prot No.:	P22676
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.29% sodium azide.
Specificity:	This antibody detects endogenous levels of Calretinin at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year