

Antibodies Product Datasheet

Product Name:	Src Mouse mAb(Mix-mA)
Catalog Number:	N1458
Reactivity:	Human, Mouse,Rat
Applications:	IHC-p,WB
Dilutions:	IHC-p1:50-200 ,WB 1:1000-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Src
Concentration:	1 mg/ml
Gene Name:	SRC SRC1
Protein Name:	Src
Other Name:	Proto-oncogene tyrosine-protein kinase Src (EC 2.7.10.2) (Proto-oncogene c-Src) (pp60c-src) (p60-Src)
Human Swiss Prot No.:	P12931
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.31% sodium azide.
Specificity:	This antibody detects endogenous levels of Src at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year