

Antibodies Product Datasheet

Product Name:	PPAR γ Mouse mAb(Mix-mA)
Catalog Number:	N1459
Reactivity:	Human, Mouse,Rat
Applications:	IHC-p,WB
Dilutions:	IHC-p1:50-200 ,WB 1:1000-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human PPAR γ
Concentration:	1 mg/ml
Gene Name:	PPARG NR1C3
Protein Name:	PPAR γ
Other Name:	Peroxisome proliferator-activated receptor gamma (PPAR-gamma) (Nuclear receptor subfamily 1 group C member 3)
Human Swiss Prot No.:	P37231
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.32% sodium azide.
Specificity:	This antibody detects endogenous levels of PPAR γ at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year