

Antibodies Product Datasheet

Product Name:	HSC70 Mouse mAb(7G12)
Catalog Number:	N1461
Reactivity:	Human, Mouse,Rat
Applications:	IHC-p,WB
Dilutions:	IHC-p1:50-200 ,WB 1:1000-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human HSC70
Concentration:	1 mg/ml
Gene Name:	HSPA8 HSC70 HSP73 HSPA10
Protein Name:	HSC70
Other Name:	Heat shock cognate 71 kDa protein (Heat shock 70 kDa protein 8)
Human Swiss Prot No.:	P11142
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.34% sodium azide.
Specificity:	This antibody detects endogenous levels of HSC70 at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year