

Antibodies Product Datasheet

Product Name:	AMPK γ 1 Mouse mAb(5H5)
Catalog Number:	N1464
Reactivity:	Human, Mouse,Rat
Applications:	IHC-p,WB
Dilutions:	IHC-p1:50-200 ,WB 1:1000-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human AMPK γ 1
Concentration:	1 mg/ml
Gene Name:	PRKAG1
Protein Name:	AMPK γ 1
Other Name:	5'-AMP-activated protein kinase subunit gamma-1 (AMPK gamma1) (AMPK subunit gamma-1) (AMPKg)
Human Swiss Prot No.:	P54619
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.37% sodium azide.
Specificity:	This antibody detects endogenous levels of AMPK γ 1 at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year