

Antibodies Product Datasheet

Product Name:	HO-1 Mouse mAb(Mix-mA)
Catalog Number:	N1466
Reactivity:	Human, Mouse,Rat
Applications:	IHC-p,WB
Dilutions:	IHC-p1:50-200 ,WB 1:1000-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human HO-1
Concentration:	1 mg/ml
Gene Name:	HMOX1 HO HO1
Protein Name:	HO-1
Other Name:	Heme oxygenase 1 (HO-1) (EC 1.14.99.3)
Human Swiss Prot No.:	P09601
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.39% sodium azide.
Specificity:	This antibody detects endogenous levels of HO-1 at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year