

Antibodies Product Datasheet

Product Name:	HIF-1 α Mouse mAb(Mix-mA)
Catalog Number:	N1473
Reactivity:	Human, Mouse,Rat
Applications:	IHC-p,WB
Dilutions:	IHC-p1:50-200 ,WB 1:1000-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human HIF-1 α
Concentration:	1 mg/ml
Gene Name:	HIF1A BHLHE78 MOP1 PASD8
Protein Name:	HIF-1 α
Other Name:	Hypoxia-inducible factor 1-alpha (HIF-1-alpha) (HIF1-alpha) (ARNT-interacting protein) (Basic-helix-loop-helix-PAS protein MOP1) (Class E basic helix-loop-helix protein 78) (bHLHe78) (Member of PAS protein 1) (PAS domain-containing protein 8)
Human Swiss Prot No.:	Q16665
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.46% sodium azide.
Specificity:	This antibody detects endogenous levels of HIF-1 α at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year