

Antibodies Product Datasheet

Product Name:	p38 Mouse mAb(4G2)
Catalog Number:	N1476
Reactivity:	Human, Mouse,Rat
Applications:	IHC-p,WB
Dilutions:	IHC-p1:50-200 ,WB 1:1000-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human p38
Concentration:	1 mg/ml
Gene Name:	MAPK14 CSBP CSBP1 CSBP2 CSPB1 MXI2 SAPK2A
Protein Name:	p38
Other Name:	Mitogen-activated protein kinase 14 (MAP kinase 14) (MAPK 14) (EC 2.7.11.24) (Cytokine suppressive anti-inflammatory drug-binding protein) (CSAID-binding protein) (CSBP) (MAP kinase MXI2) (MAX-interacting protein 2) (Mitogen-activated protein kinase p38 a
Human Swiss Prot No.:	Q16539
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.49% sodium azide.
Specificity:	This antibody detects endogenous levels of p38 at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year