

Antibodies Product Datasheet

Product Name:	Rho A Mouse mAb(4E12)
Catalog Number:	N1477
Reactivity:	Human, Mouse,Rat
Applications:	IHC-p,WB
Dilutions:	IHC-p1:50-200 ,WB 1:1000-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Rho A
Concentration:	1 mg/ml
Gene Name:	RHOA ARH12 ARHA RHO12
Protein Name:	Rho A
Other Name:	Transforming protein RhoA (Rho cDNA clone 12) (h12)
Human Swiss Prot No.:	P61586
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.50% sodium azide.
Specificity:	This antibody detects endogenous levels of Rho A at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year