

Antibodies Product Datasheet

Product Name:	Myd88 Mouse mAb(2H12)
Catalog Number:	N1479
Reactivity:	Human, Mouse,Rat
Applications:	IHC
Dilutions:	IHC 1:50-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Myd88
Concentration:	1 mg/ml
Gene Name:	MYD88
Protein Name:	Myd88
Other Name:	Myeloid differentiation primary response protein MyD88
Human Swiss Prot No.:	Q99836
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.52% sodium azide.
Specificity:	This antibody detects endogenous levels of Myd88 at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year