

Antibodies Product Datasheet

Product Name:	SOCS-1 Mouse mAb(5C3)
Catalog Number:	N1480
Reactivity:	Human, Mouse,Rat
Applications:	IHC-p,WB
Dilutions:	IHC-p1:50-200 ,WB 1:1000-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human SOCS-1
Concentration:	1 mg/ml
Gene Name:	SOCS1 SSI1 TIP3
Protein Name:	SOCS-1
Other Name:	Suppressor of cytokine signaling 1 (SOCS-1) (JAK-binding protein) (JAB) (STAT-induced STAT inhibitor 1) (SSI-1) (Tec-interacting protein 3) (TIP-3)
Human Swiss Prot No.:	O15524
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.53% sodium azide.
Specificity:	This antibody detects endogenous levels of SOCS-1 at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year