

Antibodies Product Datasheet

Product Name:	JNK2 Mouse mAb(4D2)
Catalog Number:	N1486
Reactivity:	Human, Mouse,Rat
Applications:	IHC
Dilutions:	IHC 1:50-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human JNK2
Concentration:	1 mg/ml
Gene Name:	MAPK9 JNK2 PRKM9 SAPK1A
Protein Name:	JNK2
Other Name:	Mitogen-activated protein kinase 9 (MAP kinase 9) (MAPK 9) (EC 2.7.11.24) (JNK-55) (Stress-activated protein kinase 1a) (SAPK1a) (Stress-activated protein kinase JNK2) (c-Jun N-terminal kinase 2)
Human Swiss Prot No.:	P45984
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.59% sodium azide.
Specificity:	This antibody detects endogenous levels of JNK2 at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year