

Antibodies Product Datasheet

Product Name:	NF κ B P105/P50 Mouse mAb(4D11)
Catalog Number:	N1487
Reactivity:	Human, Mouse
Applications:	IHC-p,WB
Dilutions:	IHC-p1:50-200 ,WB 1:1000-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human NF κ B P105/P50
Concentration:	1 mg/ml
Gene Name:	NFKB1
Protein Name:	NF κ B P105/P50
Other Name:	Nuclear factor NF-kappa-B p105 subunit (DNA-binding factor KBF1) (EBP-1) (Nuclear factor of kappa light polypeptide gene enhancer in B-cells 1) [Cleaved into: Nuclear factor NF-kappa-B p50 subunit]
Human Swiss Prot No.:	P19838
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.60% sodium azide.
Specificity:	This antibody detects endogenous levels of NF κ B P105/P50 at Human, Mouse
Storage and Stability:	-20°C/1 year

N1487 | NF κ B P105/P50 Mouse mAb(4D11)

Tech@assaybiotech.com | Order@assaybiotech.com

Research use only not intended for therapeutic/diagnostic procedures

Last updated December 2022