

Antibodies Product Datasheet

Product Name:	Yeast GAPDH Mouse mAb(Mix-mA)
Catalog Number:	N1493
Reactivity:	Yeast (H, M, R)
Applications:	IHC-p,WB
Dilutions:	IHC-p1:50-200 ,WB 1:1000-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from Yeast GAPDH
Concentration:	1 mg/ml
Gene Name:	Yeast GAPDH
Protein Name:	Yeast GAPDH
Other Name:	
Human Swiss Prot No.:	P60010
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.66% sodium azide.
Specificity:	This antibody detects endogenous levels of Yeast GAPDH
Storage and Stability:	-20°C/1 year