

Antibodies Product Datasheet

Product Name:	Drosophila Melanogaster GAPDH Mouse mAb(7B7)
Catalog Number:	N1496
Reactivity:	Drosophila Melanogaster (H, M, R)
Applications:	IHC-p, WB
Dilutions:	IHC-p1:100-500, WB 1:5000-20000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from Drosophila Melanogaster GAPDH
Concentration:	1 mg/ml
Gene Name:	Drosophila Melanogaster GAPDH
Protein Name:	Drosophila Melanogaster GAPDH
Other Name:	
Human Swiss Prot No.:	P07486
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.69% sodium azide.
Specificity:	This antibody detects endogenous levels of Drosophila Melanogaster GAPDH at Drosophila Melanogaster (H, M, R)
Storage and Stability:	-20°C/1 year