

Antibodies Product Datasheet

Product Name:	HSP60 Mouse mAb(2D3)
Catalog Number:	N1500
Reactivity:	Human, Mouse,Rat
Applications:	IHC-p,WB
Dilutions:	IHC-p1:50-200 ,WB 1:1000-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human HSP60
Concentration:	1 mg/ml
Gene Name:	HSPD1 HSP60
Protein Name:	HSP60
Other Name:	60 kDa heat shock protein, mitochondrial (60 kDa chaperonin) (Chaperonin 60) (CPN60) (Heat shock protein 60) (HSP-60) (Hsp60) (HuCHA60) (Mitochondrial matrix protein P1) (P60 lymphocyte protein)
Human Swiss Prot No.:	P10809
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.73% sodium azide.
Specificity:	This antibody detects endogenous levels of HSP60 at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year