

Antibodies Product Datasheet

Product Name:	JAK1 Mouse mAb(3F12)
Catalog Number:	N1502
Reactivity:	Human, Mouse,Rat
Applications:	IHC-p,WB
Dilutions:	IHC-p1:50-200 ,WB 1:1000-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human JAK1
Concentration:	1 mg/ml
Gene Name:	JAK1 JAK1A JAK1B
Protein Name:	JAK1
Other Name:	Tyrosine-protein kinase JAK1 (EC 2.7.10.2) (Janus kinase 1) (JAK-1)
Human Swiss Prot No.:	P23458
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.78% sodium azide.
Specificity:	This antibody detects endogenous levels of JAK1 at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year

N1502 | JAK1 Mouse mAb(3F12)

Tech@assaybiotech.com | Order@assaybiotech.com

Research use only not intended for therapeutic/diagnostic procedures

Last updated December 2022